

deconex® 12 PA-x

High performance cleaning concentrate for use in research, medicine, biology and chemistry as well as for optical and electronic components

For soaking baths and ultrasonic cleaning systems



Application

deconex 12 PA-x is a liquid, mildly alkaline yet highly efficient concentrate for the spontaneous cleaning and decontamination of delicate laboratory equipment and as well of optical and electronic components (e. g. semiconductors) out of glass, quartz, porcelain, metal etc.

Properties

deconex 12 PA-x removes the most stubborn dirt such as burnt-in grease, oil, wax, dried-on blood, proteins, residues of cell cultures or paint, pigments, organic and inorganic traces as well as silicone.

Dosage

The optimum dosage depends mainly on the hardness of the water and the stubbornness of the contamination. Using demineralized water improves the cleaning properties.

In practice, the following dosages have given good results:

Procedure	Dosage	Temperature	Exposure time
Normal	0.5-1.0 %	20-30 °C	1-2 h
Intensive	2.0-4.0 %	20-30 °C	6-12 h
Fast*	2.0-5.0 %	40-70 °C	2-5 min

* with ultrasonic cleaner

The cleaning procedure is accelerated by raising the temperature of the solution. The time necessary is roughly halved for each increase of 10 °C.

deconex 12 PA-x considerably reduces the surface tension of the cleaning solution and, because of its wetting performance, removes dirt easily. The excellent emulsifying and dispersing capabilities prevent removed particles from being redeposited. Special surfactants guarantee excellent rinsability of the cleaned surface. Instead of etching and corrosive substances such as caustic alkalis and chlorine, mild cleaning enhancers have been used to make the product safe for materials.

deconex 12 PA-x is NTA-free.

Ingredients

Alkalis, complexing agents, dispersing agent, solubiliser, surface active compounds, wetting agent

deconex® 12 PA-x

Material compatibility

Suitable for:
Lead crystal, steel, glass, rubber, ceramic, porcelain, quartz,
chrome steel, PE, PP, PVC, PTFE, Viton, NBR

For materials not mentioned please make your own specific compatibility tests or consult Borer Chemie AG.

Chemical/physical data

pH	1 %	approx. 11.7
Density	concentrate	1.42 g/ml
Appearance	concentrate	light yellowish

Availability

Please ask your local representative about current container sizes.

Containers, screw caps and labels are made of recyclable polyethylene.

Additional information

For information concerning safety at work, storage and waste disposal/effluent, please consult the corresponding safety data sheet.

Take advantage of our vast know-how! Please, contact us for further information regarding your specific application.

Manufacturer:

Borer Chemie AG

Gewerbstrasse 13, 4528 Zuchwil / Switzerland

Tel +41 32 686 56 00 Fax +41 32 686 56 90

office@borer.ch, www.borer.ch

All information provided is based on our current knowledge and it does not constitute a legally binding assurance of specific product properties.