

## deconex® dosing table for disinfectants and detergent concentrates

Instructions to prepare a ready-to-use solution of a deconex® disinfectant or detergent. The quantity of use solution consists of tap water and the respective disinfectant or detergent concentrate (e.g. to receive 10 litres of a 2% solution, 200 ml of concentrate will have to be added to 9.8 litres of tap water). Follow manufacturer's recommendations and guidelines.

Percentage of ready-to-use solution	0.25 %	0.50 %	0.75 %	1.00 %	1.50 %	2.00 %	3.00 %	4.00 %	5.00 %
Total volume of ready-to-use solution	Required amount of disinfectant or detergent concentrate (ml)								
<b>0.5 litre</b>	1.25	2.5	3.75	5	7.5	10	15	20	25
<b>1 litre</b>	2.5	5	7.5	10	15	20	30	40	50
<b>2 litre</b>	5	10	15	20	30	40	60	80	100
<b>3 litre</b>	7.5	15	22.5	30	45	60	90	120	150
<b>4 litre</b>	10	20	30	40	60	80	120	160	200
<b>5 litre</b>	12.5	25	37.5	50	75	100	150	200	250
<b>6 litre</b>	15	30	45	60	90	120	180	240	300
<b>7 litre</b>	17.5	35	52.5	70	105	140	210	280	350
<b>8 litre</b>	20	40	60	80	120	160	240	320	400
<b>9 litre</b>	22.5	45	67.5	90	135	180	270	360	450
<b>10 litre</b>	25	50	75	100	150	200	300	400	500
<b>15 litre</b>	37.5	75	112.5	150	225	300	450	600	750
<b>20 litre</b>	50	100	150	200	300	400	600	800	1000
<b>25 litre</b>	62.5	125	187.5	250	375	500	750	1000	1250
<b>30 litre</b>	75	150	225	300	450	600	900	1200	1500