

deconex® 33 COMPLEX

Liquid complexing agent

For use in laboratory glassware washers



Application

deconex 33 COMPLEX is a liquid, complexing additive for use in the presence of metallic salt residues or heavy metal ions or with especially hard water.

deconex 33 COMPLEX can be used as an additive to alkaline cleaners in the main washing cycle or as an additive to the rinsing water.

Properties

deconex 33 COMPLEX is phosphate-free and has a good complexing capacity for metal ions, especially for alkaline earth and heavy metal ions.

Ingredients

Complexing agent

Dosage

As an additive to alkaline cleaners or in rinsing water

0.2-1 ml/l

Chemical/physical data

pH	1%	approx. 11.0
Density	concentrate	1.30 g/ml
Appearance	concentrate	colourless

Availability

Please ask your local representative about current container sizes.

Containers, screw caps and labels are made of recyclable polyethylene.

Additional information

For information concerning safety at work, storage and waste disposal/effluent, please consult the corresponding safety data sheet.

Take advantage of our vast know-how! Please, contact us for further information regarding your specific application.

deconex® 33 COMPLEX

Manufacturer:

Borer Chemie AG

Gewerbestrasse 13, 4528 Zuchwil / Switzerland

Tel +41 32 686 56 00 Fax +41 32 686 56 90

office@borer.ch, www.borer.ch

All information provided is based on our current knowledge and it does not constitute a legally binding assurance of specific product properties.



advanced cleaning solutions