



deconex® 41 STS ALKA

The new triple alkaline cleaner

The deconex® 41 STS product family has been enlarged with a new product that offers 3 solutions in one! The new deconex® 41 STS ALKA is an alkaline agent that prevents lime scale and acts as a cleaning and rinsing agent in automated decontamination machines for human waste containers.

deconex® 41 STS ALKA has particularly been developed for machines containing acid-sensitive components and/or a dosing system which injects the product directly into the steam generator. Its complexing agent prevents formation of boiler scale in the steam generator. The mild alkaline product also prevents formation of deposits on accessories, machine walls or in nozzles and pipes.

The new deconex® 41 STS ALKA offers the following advantages

Efficient

- Excellent cleaning, rinsing and descaling power
- Prevents formation of deposits in steam generator and on washing goods
- Suitable for stainless steel, copper, cast iron, plastics (e.g. polycarbonate bottles)

Safe

- Certified quality product
- Compliance with highest hygienic requirements
- Extensively tested efficacy, safety and material compatibility
- Environmental-friendly, biodegradable product

Economical

- Cost effective - one product for 3 applications
- Synergetic effect of quality ingredients for highest efficiency at minimum dosage
- Suitable for all water qualities

The deconex® 41 STS product family



The deconex® 41 STS product family overview

The deconex® 41 STS products have especially been developed for the cleaning of equipment like bed pans, urine bottles and other human waste containers.

Cleaning	Rinsing	Descaling
deconex® 41 STS ALKA	deconex® 41 STS ALKA	deconex® 41 STS ALKA
	deconex® 41 STS	deconex® 41 STS
	deconex® 41 STS TF	deconex® 41 STS TF

Please visit our website or contact us or your local distributor for further information.

Borer Chemie AG

Gewerbstrasse 13, 4528 Zuchwil / Switzerland
 Tel. +41 32 686 56 00, Fax +41 32 686 56 90
 www.borer.ch, office@borer.ch

